

# Specifications for the ENOCHS EXCEL 355 General Examination Table

<b>Warranty:</b>	4-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Back Lift Mechanism:</b>	Horizontal to 80° by spring assisted counter balance
<b>Couch Top:</b>	One-piece upholstered top, steel reinforced
<b>Perineal Recess:</b>	Not available on this model
<b>Pelvic Tilt:</b>	Not available on this model
<b>Leg Extension:</b>	Extends out 18" from end of seat section
<b>Irrigation Basin:</b>	4" Deep basin
<b>Electrical Outlets:</b>	Dual electrical outlets, hospital grade
<b>Instrument Warmer: with ON/OFF Switch, Circuit Breaker &amp; Reset</b>	Not available on this model
<b>Stirrups:</b>	Multi-positional, lateral and length adjustment
<b>Storage:</b>	<b>3 Side Drawers:</b> 21.5" W x 15" D x 4" H <b>2 Footend Drawers:</b> 15.5" W x 15" D x 4" H <b>Total Storage Space:</b> 3.32 cubic feet
<b>Foot Step:</b>	Extra long surface engineered for 600 lb. patient load
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	knee crutches, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with drain and screen, footend accessory rails, caster base
<b>Mains Power Entry: with Circuit Breaker</b>	Standard feature
<b>Duty Cycle:</b>	<b>Receptacles:</b> Continuous <b>Instrument Warming Drawer:</b> n/a
<b>Ratings:</b>	<b>Table:</b> 115V 60Hz 10 amps <b>Receptacles:</b> 115V 60Hz <b>Total Receptacle Load:</b> 9 amps
<b>Table Dimensions:</b>	<b>Width:</b> 27", <b>Height:</b> 34", <b>Length:</b> 55" (extends to 73")
<b>Table Weight:</b>	268 lbs.
<b>Patient Load:</b>	600 lbs. evenly distributed on main couch
<b>Certifications:</b>	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1