

# Specifications for the RESOLUTE 6220

<b>Warranty:</b>	4-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Power System:</b>	Electro-mechanical
<b>Attainable Positions:</b>	Chair, Contoured Chair, or OB/GYN, Lithotomy, Urology, Flat
<b>Power Adjustments:</b>	<b>Height:</b> 24" to 40", <b>Back:</b> horizontal 80° <b>Tilt:</b> horizontal to 30°, <b>Leg:</b> vertical to horizontal
<b>Programmable:</b>	Programmable with optional hand control accessory
<b>Hand Control:</b>	Optional accessory that may be added later
<b>Hand Control Connection Port:</b>	Standard feature, located on both sides of table's moveable column for easy access and elimination of over-stretched cords.
<b>Foot Control:</b>	"Easy touch" hands-free foot control Safety Extra-Low Voltage (SELV) and IPX1
<b>Upholstered Couch Top:</b>	One-piece style, no cross seams, embossing or folded fabric, securely fastened to all steel frame
<b>Headend Accessory Rails:</b>	Stainless steel rails, mounted on both sides of the table
<b>Irrigation Basin:</b>	2" deep stainless steel basin, integrated safety switch prohibits operation of foot section when basin is extended
<b>Electrical Receptacles:</b>	Four hospital grade electrical receptacles (2 per side) Overload Protection: Dedicated circuit breaker with quick reset
<b>Safety Lockout:</b>	Shuts power off to all articulating sections of table, but retains power to receptacles
<b>Stirrups:</b>	Positive lock, multi-position stirrups, length and lateral adjustment
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	Knee crutches, padded arm board, restraint straps, IV pole, urology basin with screen and drain, footend accessory rails, caster base, universal clamp, multi-positional clamp
<b>Mains Power Entry:</b>	Master ON/OFF, IEC connection, Overload Protection: Integrated circuit breaker with quick reset
<b>Duty Cycle:</b>	1 minute ON / 9 minutes OFF
<b>Ratings:</b>	<b>Table:</b> 115V 60 Hz 10 amps <b>Receptacles:</b> 115V 60 Hz <b>Total Receptacle Load:</b> 3.9 amps
<b>Table Dimensions:</b>	<b>Length:</b> 71", extends with headrest to 83" <b>Width top:</b> 30" at seat, tapers to 25" at headend <b>Headrest:</b> 19" wide, tapers to 15" <b>Legrest:</b> 19" wide, tapers to 18"
<b>Table Weight:</b>	485 lbs.
<b>Patient Load:</b>	500 lbs. evenly distributed on main couch
<b>Certifications:</b>	UL-60601-1, Can/CSA C22.3 No. 601-1