

# ENOCHS ENCORE 4200

## High-Low Examination Table

**PATIENT LOAD: 500 LBS.**  
evenly distributed on the table's main couch

**PAPER ROLL HOLDER** is concealed at headend of table and accommodates 18" or 21" paper

**PAPER CUTTER** positioned at the headend of table, encourages staff to remove the entire length of table paper contaminated by patient

**FULL POWER CONTROLLED BACK** smoothly adjusts to any angle from horizontal to an 80° chair position

**FIELD REVERSIBLE SIDE DRAWER**

Durable 6" deep **MOULDED LARGE SIDE DRAWER** with rounded internal corners and edges for easy and efficient cleaning and maintenance

Extra strength **STEEL REINFORCED SIDE GUIDE MOUNTING SYSTEM** for discipline non-binding movement

The 100% seamless **VACUUM-FORMED TOP** is designed without decorative embossing, cross-seams and fabric folds for easy and efficient cleaning. Top is padded with 1.5" of high-density foam with an extra .5" along the edges for patient comfort.

**PELVIC TILT** - 7° and 11° positions with release handles on either side of table

100% **PERINEAL RECESS** seamless for infection control maintenance

4" deep **STAINLESS STEEL IRRIGATION BASIN** moves independent of leg extension

**POSITIVE LOCK MULTI-POSITION STIRRUPS** with length and lateral adjustment

**PADDED PULLED LEG EXTENSION** works independent of Irrigation Basin and allows full access to Basin when unit is extended

4" Deep **DRAWER** with rounded internal corners and edges for quick and easy cleaning

**MAINS POWER LIGHT**

**'EASY TOUCH' FOOT CONTROL** adjusts table's height and back with reliable hands-free action. Pedals are covered with non-slip tread and utilize a natural left to right rocking motion of operator's foot.

### TABLE FEATURES:

- 2-YEAR Warranty
- Hand Control available as an accessory
- Fully welded All Steel Base Construction - NO PLASTIC
- Full Power adjustment of back
- Easy to clean 100% Seamless Vacuum Formed Top with additional padding on edges for patient comfort
- Pelvic Tilt 7° and 11°
- Stainless Steel Irrigation Basin operates independent of Leg Extension
- Easy to adjust Stirrups with length and lateral adjustment
- 1.4 cubic feet of Storage
- Certified by UL to Standard 60601-1 CAN/CSA C22.3 No. 601-1
- ENOCHS Encore 4200 replaces the Power 4000-01



Securely attached to the table base, the ENOCHS vacuum-formed top is the strongest top on the market with an all steel frame and back. Top is padded for patient comfort with 1.5" of high density foam with an additional .5" of foam contouring the edges. The sleek Euro-look top features a seamless surface, moulded perineal recess and a pelvic tilt with positions of 7° and 11°. The top offers proven performance of strength, durability and long service life.



The working end of the table has the exclusive ENOCHS design of the pullout padded leg extension moving independent of the 4" deep stainless steel irrigation basin. This design reduces surface area to allow fast easy clean up after irrigation procedure. Patients with limited rotation can be placed in a GYN position utilizing the pad board instead of stirrups while physician maintains unrestricted access to the irrigation basin.

# Specifications for the ENOCHS ENCORE 4200 High-Low Examination Table

<b>Warranty:</b>	2-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Power System:</b>	Electro-Mechanical provides smooth and quiet movement
<b>Power Adjustments:</b>	<b>Height:</b> 24" - 40" <b>Back:</b> horizontal - 80 °
<b>Foot Control:</b>	"Hands free", standard on every table, SELV
<b>Hand Control:</b>	"Easy to read", optional accessory, add at anytime, SELV
<b>Couch Top:</b>	100% Seamless vacuum formed style, all steel frame construction
<b>Perineal Recess:</b>	Moulded, no seams, 10" W x 2" D at center point
<b>Pelvic Tilt:</b>	7°, 11°
<b>Leg Extension:</b>	Extends out 20" from end of seat section
<b>Irrigation Basin:</b>	4" Deep stainless steel basin
<b>Electrical Outlets: with Circuit Breaker</b>	Not available on this model
<b>Instrument Warmer: with ON/OFF Switch, Circuit Breaker &amp; Reset</b>	Not available on this model
<b>Stirrups:</b>	Multi-positional, lateral and length adjustment
<b>Storage:</b>	<b>Side Drawer:</b> 15.5" W x 15" D x 6" H <b>Footend Drawer:</b> 15.5" W x 15" D x 4" H <b>Total Storage Space:</b> 1.4 cubic feet
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	hand control, knee crutches, padded arm board, drawer locks, restraint straps, IV pole, urology basin with screen and drain, headend accessory rails, footend accessory rails, caster base
<b>Mains Power Entry: with Circuit Breaker</b>	IEC connection, internal circuit breaker with reset
<b>Duty Cycle:</b>	<b>Table:</b> 1 minute ON / 9 minutes OFF <b>Receptacles:</b> n/a <b>Instrument Warming Drawer:</b> n/a
<b>Ratings:</b>	<b>Table:</b> 115V 60 Hz 10 amps <b>Receptacles:</b> n/a <b>Total Receptacle Load:</b> n/a
<b>Table Dimensions:</b>	<b>Width:</b> 27", <b>Height:</b> 24" - 40", <b>Length:</b> 52" (extends to 72")
<b>Table Weight:</b>	416 lbs.
<b>Patient Load:</b>	500 lbs. evenly distributed on main couch
<b>Certifications:</b>	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1