

ENOCHS EXCEL 355

General Examination Table

PATIENT LOAD: 600 LBS.
evenly distributed on the table's main couch

SPRING ASSISTED COUNTER BALANCE MECHANISM for smooth easy adjustment of the back section - NO gas cylinder to replace

ONE-PIECE UPHOLSTERED TOP padded with 1.5" high density foam on a steel reinforced frame for maximum strength

Concealed when not in use, the **MULTI-POSITION STIRRUPS** offer length and lateral adjustment

PAPER ROLL HOLDER

PAPER CUTTER

MAINS POWER LIGHT illuminates when table has power

EXTRA WIDE PULL-OUT LEG EXTENSION

FIELD REVERSIBLE SIDE DRAWERS

4" deep IRRIGATION BASIN moves independent of leg extension

Durable **4" DEEP MOULDED DRAWERS** with rounded edges and corners for easy cleaning

FOOTEND STORAGE DRAWERS

Easy to clean **SLIP-RESISTANT TREAD**

Extra strength **STEEL REINFORCED SIDE GUIDE MOUNTING SYSTEM** for discipline non-binding movement

PULL-OUT LOCKING STEP engineered to accommodate the 600 lbs patient load

FULLY WELDED STEEL BASE construction – NO PLASTIC PARTS

TABLE FEATURES:

- 4-YEAR Warranty
- All Steel Base Construction - NO PLASTIC
- Dual Electrical Outlets, Hospital Grade, 9 amps total load
- Easy to clean One-Piece Upholstered Top, designed without cross-seams, decorative stitching or embossing
- Top constructed with a Steel Reinforced Frame for maximum strength
- Irrigation Basin moves independent of Leg Extension for fast transition
- Multi-position Stirrups with easy length and lateral adjustment
- 3.32 cubic feet of Storage in the five spacious Drawers
- Certified by UL to Standard 60601-1 CAN/CSA C22.3 No. 601-1
- ENOCHS Excel 355 replaces the Excel 71



The extra long pull-out foot step features a steel stabilizing leg designed to accommodate a 600 pound patient load capacity. The step easily extends out from the table and locks into position providing the patient with a stable surface when getting on and off the table. The step's all steel construction offers extraordinary strength and durability. The 17.5" long surface is covered in easy to clean slip-resistant tread.



New field reversible side drawers allows table to have the flexibility to be placed into any room configuration. Side drawers combined with the two footend drawers offer maximum storage capacity totaling 3.32 cubic feet. Nylatron® on the bottom drawer edges protect against wear and enhances the drawer's overall performance. An integrated stops eliminate accidental spillage of drawer contents.



The table's electrical receptacles feature a new tool required access door designed to restrict unauthorized use of the outlets. Door easily unlocks so accessories can be plugged into the UL Listed, hospital grade receptacles. Once in place, the door closes over the electrical cords, securing plugs in place with the locking mechanism. Receptacles have a 9 amps total maximum load.

Specifications for the ENOCHS EXCEL 355 General Examination Table

Warranty:	4-YEAR
Base Construction:	Fully welded all steel base
Back Lift Mechanism:	Horizontal to 80° by spring assisted counter balance
Couch Top:	One-piece upholstered top, steel reinforced
Perineal Recess:	Not available on this model
Pelvic Tilt:	Not available on this model
Leg Extension:	Extends out 18" from end of seat section
Irrigation Basin:	4" Deep basin
Electrical Outlets:	Dual electrical outlets, hospital grade
Instrument Warmer: with ON/OFF Switch, Circuit Breaker & Reset	Not available on this model
Stirrups:	Multi-positional, lateral and length adjustment
Storage:	3 Side Drawers: 21.5" W x 15" D x 4" H 2 Footend Drawers: 15.5" W x 15" D x 4" H Total Storage Space: 3.32 cubic feet
Foot Step:	Extra long surface engineered for 600 lb. patient load
Paper Roll Holder:	Accommodates 18" or 21" paper widths
Accessories:	knee crutches, drawer locks, stainless steel basin, restraint straps, IV pole, urology basin with drain and screen, footend accessory rails, caster base
Mains Power Entry: with Circuit Breaker	Standard feature
Duty Cycle:	Receptacles: Continuous Instrument Warming Drawer: n/a
Ratings:	Table: 115V 60Hz 10 amps Receptacles: 115V 60Hz Total Receptacle Load: 9 amps
Table Dimensions:	Width: 27", Height: 34", Length: 55" (extends to 73")
Table Weight:	268 lbs.
Patient Load:	600 lbs. evenly distributed on main couch
Certifications:	By UL to Standard 60601-1, CAN/CSA C22.3 No. 601-1