

# ENOCHS RESOLUTE 6210

## Basic Procedural Table

**PATIENT LOAD: 500 LBS.**  
evenly distributed on the table's main couch



### TABLE FEATURES:

- **2-YEAR WARRANTY**
- "Easy Touch" Foot Control for hands free adjustment of table
- One-piece Upholstered Top designed without cross seams, decorative embossing or material folds that hinder inflectional control procedures
- 30" Wide Seat Section to provide larger patients more comfort
- Mains ON/OFF shuts off all power to table to prevent accidental movement of table during exam or procedure
- Mains Power Overload Protection: easy reset Circuit Breaker - no fuses
- Adjustable Padded Headrest
- Multi-Position Positive Lock Stirrups
- 2" deep Stainless Steel Irrigation Basin
- Stainless Steel Headend Accessory Rails



Positive lock, multi-position stirrups with lateral and length adjustment and a 2" deep stainless steel irrigation basin that quickly and easily removes for cleaning after each procedure are standard features of the 6210. Table is engineered with an integrated safety switch that stops the leg section from going UP when irrigation basin is not in full stored position, eliminating damage to the basin and/or leg extension.



The table's 30" wide seat section provides maximum comfort to larger patients while the one-piece seamless top surface, designed without decorative embossing or material creases and folds that collect dirt, allows staff to quickly and efficiently clean after each patient. The fully padded upholstered top, covered in durable Naugahyde®, is securely attached to an all steel frame for strength and stability.

# Specifications for the RESOLUTE 6210

<b>Warranty:</b>	2-YEAR
<b>Base Construction:</b>	Fully welded all steel base
<b>Power System:</b>	Electro-mechanical
<b>Attainable Positions:</b>	Chair, Contoured Chair, Trendelenburg, OB/GYN, Lithotomy, Urology, Flat
<b>Power Adjustments:</b>	<b>Height:</b> 24" to 40", <b>Back:</b> horizontal to 80° <b>Tilt:</b> horizontal to 30°, <b>Leg:</b> vertical to horizontal
<b>Programmable:</b>	Not available on this model
<b>Hand Control:</b>	Not available on this model
<b>Hand Control Connection Port:</b>	Not available on this model
<b>Foot Control:</b>	"Easy touch" hands-free foot control Safety Extra-Low Voltage (SELV) and IPX1
<b>Upholstered Couch Top:</b>	One-piece style, no cross seams, embossing or folded fabric, securely fastened to all steel frame
<b>Headend Accessory Rails:</b>	Stainless steel rails, mounted on both sides of the table
<b>Irrigation Basin:</b>	2" deep stainless steel basin, integrated safety switch prohibits operation of foot section when basin is extended
<b>Electrical Receptacles:</b>	Not available on this model
<b>Safety Lockout:</b>	Not available on this model
<b>Stirrups:</b>	Positive lock, multi-position stirrups, length and lateral adjustment
<b>Paper Roll Holder:</b>	Accommodates 18" or 21" paper widths
<b>Accessories:</b>	Knee crutches, padded armboard, restraint straps, IV pole, urology basin with screen and drain, footend accessory rails, caster base, universal clamp, multi-positional clamp
<b>Mains Power Entry:</b>	Master ON/OFF, IEC connection, Overload Protection: Integrated circuit breaker with quick reset
<b>Duty Cycle:</b>	1 minute ON / 9 minutes OFF
<b>Ratings:</b>	<b>Table:</b> 115V 60 Hz 10 amps <b>Receptacles:</b> N/A, <b>Total receptacle load:</b> N/A
<b>Table Dimensions:</b>	<b>Length:</b> 71", extends with headrest to 83" <b>Width top:</b> 30" at seat, tapers to 25" at headend <b>Headrest:</b> 19" wide, tapers to 15" <b>Legrest:</b> 19" wide, tapers to 18"
<b>Table Weight:</b>	485 lbs.
<b>Patient Load:</b>	500 lbs. evenly distributed on main couch
<b>Certifications:</b>	UL 60601-1, Can/CSA C22.3 No. 601-1